

CONVENIENCE
COMFORT
HYGIENIC PROTECTION

DENTAL CARE

Single-use infection
prevention products

 **TIDI**
PRODUCTS
Care with Confidence



Care with Confidence

TIDI focuses on you, the healthcare professional, continually evaluating your needs. So we understand that it's imperative to create an exam environment that presents a professional image and makes patients feel comfortable. It's essential that exam products be easy and convenient to use. And it's critical to minimize any chance for cross contamination.

In response, we engineer specific single-use infection prevention solutions to meet these needs. You can be confident in our reliable, consistent product quality, assured by our certification in ISO 9001 and ISO 13485.

With these assurances, healthcare professionals confidently reach for TIDI's single-use infection control products to maintain a hygienic environment, put staff and patients at ease, and focus on care.



Clear Choices for Confident Care

Our simplified product organization includes three distinct levels that clearly describe our differentiated, value-added products, based on insight gleaned from healthcare professionals. This makes it easy to choose the best TIDI products for your procedures, your environment, and the comfort and convenience of you and your patients.

PRODUCT LEVELS

U Ultimate

High-performance products designed with thoughtful features for the highest-quality care, making the professional's job easier while enhancing the patient experience.

C Choice

A wide variety of convenient, high-quality products that offer comfort and protection on demand.

E Everyday

Economical single-use products that provide functional everyday protection.

TABLE OF CONTENTS

2-7	Bibs + Towels
8-9	Headrest Covers
10-11	Tray Covers
12-15	Sheaths + Barriers
16-19	Eyewear + Facemasks
20-21	Cups + Saliva Ejectors
22	Cotton + Gauze + Nonwovens
23	Sterilization + Autoclave
24	Disposal + Office Sundry
25	Glossary

Sustainable Solutions

TIDI Products is focused on reducing resource consumption, minimizing environmental impact and acting in an environmentally responsible manner, in order to provide effective solutions without compromising the well-being of future generations.

Product Samples

Experience the softness of a headrest cover, the absorbency of a towel, the convenience of a barrier, or the performance of any TIDI product with a no-charge sample pack. **Call 1.800.521.1314 or visit tidiproductions.com.**

Online Information

Explore our broad range of product areas and learn more about how our solutions can enhance the care experience. Visit tidiproductions.com.


Value-Added Packaging

Newly designed packaging provides clear, concise, instant communication. Quickly and simply verify product identity and description for ease in ordering and selecting the right TIDI products for your needs.



BIBS + TOWELS

We understand that different procedures require different kinds of protection for your patients. That's why we offer bibs with varied absorbencies and sizes to meet your needs.



TIDI Ultimate 917400
Lavendar 3-ply tissue,
poly-back, waffle-
embossed bib


U Ultimate

BIBS + TOWELS


Our most absorbent choices are excellent for prophy and minor procedures when you need to protect patients from fluids.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
919401	White	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919402	Green	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919403	Blue	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919404	Yellow	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919405	Gray	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919406	Mauve	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919407	Beige	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919408	Peach	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919410	Teal	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
919400	Lavender	13" x 18"	500	3-ply tissue, poly-back, diamond-embossed
917401	White	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917402	Green	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917403	Blue	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917404	Yellow	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917405	Gray	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917406	Mauve	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917407	Beige	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917408	Peach	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917410	Teal	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917400	Lavender	13" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917411	White	17" x 18"	500	3-ply tissue, poly-back, waffle-embossed
917413	Blue	17" x 18"	500	3-ply tissue, poly-back, waffle-embossed

Ultimate COLORS



Ultimate TEXTURES



C Choice

BIBS + TOWELS

The right bibs to choose for reliable patient protection, suited for exams, consultations and x-rays.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
919461	White	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919462	Green	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919463	Blue	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919464	Yellow	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919465	Gray	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919466	Mauve	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919467	Beige	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919468	Peach	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919470	Teal	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
919459	Lavender	13" x 18"	500	2-ply tissue, poly-back, diamond-embossed
917461	White	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917462	Green	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917463	Blue	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917464	Yellow	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917465	Gray	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917466	Mauve	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917467	Beige	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917468	Peach	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917470	Teal	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917459	Lavender	13" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917471	White	17" x 18"	500	2-ply tissue, poly-back, waffle-embossed
917473	Blue	17" x 18"	500	2-ply tissue, poly-back, waffle-embossed

Choice COLORS



Choice TEXTURES



C Choice BIBS + TOWELS

TIDI offers a variety of relaxing and soothing prints for patients of all ages.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
983813	Cat in the Hat [†]	13" x 18"	500	Specialty printed 2-ply tissue, poly-back, waffle-embossed
981410	Bugs & Things	10" x 13"	250	Specialty printed 2-ply tissue, poly-back, waffle-embossed
981413	Bugs & Things	13" x 18"	250	Specialty printed 2-ply tissue, poly-back, waffle-embossed
981913	Schooltime	13" x 18"	500	Specialty printed 2-ply tissue, poly-back, waffle-embossed
917490	TIDI Tooth	13" x 18"	500	Specialty printed 2-ply tissue, poly-back, waffle-embossed
981610	Under the Sea	10" x 13"	250	Specialty printed 2-ply tissue, poly-back, waffle-embossed
981613	Under the Sea	13" x 18"	500	Specialty printed 2-ply tissue, poly-back, waffle-embossed
982413	Watercolors	13" x 18"	500	Specialty printed 2-ply tissue, poly-back, waffle-embossed
982433	Watercolors	13" x 18"	500	Specialty printed 3-ply tissue, poly-back, waffle-embossed
983910	Horton Hears a Who [†]	10" x 13"	250	Specialty printed 2-ply tissue, poly-back, waffle-embossed
983913	Horton Hears a Who [†]	13" x 18"	500	Specialty printed 2-ply tissue, poly-back, waffle-embossed

[†]Dr. Seuss Properties TM & © 2004 Dr. Seuss Enterprises, L.P. All Rights Reserved.

Choice PATTERNS



Watercolor

Bugs & Things

Cat in the Hat[†]Horton Hears a Who[†]

Schooltime

Under the Sea

TIDI Tooth

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
7400	Clear	7" x 5" x 6"	1	Acrylic countertop dispenser fits 13" x 18" bibs and towels



Convenient acrylic countertop dispenser fits all TIDI 13" x 18" bibs and towels.



SPECIALTY BIBS



U Ultimate SPECIALTY BIBS

The highest level of absorbency for procedures such as root canals, extractions, and ultrasonic scaling. The contour neck assures a perfect fit and comfort for patients.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
917901	White	18" x 22"	400	Contour bib, 3-ply tissue, poly-back
917902	Green	18" x 22"	400	Contour bib, 3-ply tissue, poly-back
917903	Blue	18" x 22"	400	Contour bib, 3-ply tissue, poly-back
917906	Mauve	18" x 22"	400	Contour bib, 3-ply tissue, poly-back
917909	Lavender	18" x 22"	400	Contour bib, 3-ply tissue, poly-back
917441	White	17" x 18"	500	Contour bib, 3-ply tissue, poly-back
917443	Blue	17" x 18"	500	Contour bib, 3-ply tissue, poly-back
917677	Blue	18" x 25"	250	Oral surgery bib, 3-ply tissue, poly-back, shell pattern
917678	Mauve	18" x 25"	250	Oral surgery bib, 3-ply tissue, poly-back, shell pattern

C Choice SPECIALTY BIBS

High absorbency plus a contour neck provide patients with confidence and comfort during exams.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
917341	White	17" x 18"	500	Contour bib, 2-ply tissue, poly-back
917343	Blue	17" x 18"	500	Contour bib, 2-ply tissue, poly-back
917344	Mauve	17" x 18"	500	Contour bib, 2-ply tissue, poly-back
917345	Lavender	17" x 18"	500	Contour bib, 2-ply tissue, poly-back



COLORS



White



Green



Blue



Mauve




Lavender

E Everyday SPECIALTY BIBS

A selection of bibs providing substantial coverage and poly barriers to help protect patients from bleach and fluids during endodontic and oral surgery procedures.

CODE	COLOR/ PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
917673	Blue	18" x 25"	250	Chainless bib with ties, tissue/poly
917676	Mauve	18" x 25"	250	Chainless bib with ties, tissue/poly
917480	Lavender	29" x 21"	100	Endo bib, waist-length with contour neck, tissue/poly/tissue
917483	Blue	29" x 21"	100	Endo bib, waist-length with contour neck, tissue/poly/tissue
917486	Mauve	29" x 21"	100	Endo bib, waist-length with contour neck, tissue/poly/tissue
917493	Blue	29" x 42"	50	Oral surgery bib, knee-length with contour neck, tissue/ poly/tissue
917670	White	36" x 40"	100	Patient apron, cutout neck, poly

Everyday COLORS



917676 Tissue/poly bib with ties



Towels and bibs in a wide range of colors to enhance any office décor.



HEADREST COVERS

Protect your patients and equipment with our headrest covers. Each selection offers distinct characteristics to meet the needs of your office environment.



TIDI Choice 919811
White 13" x 10"
tissue/poly headrest cover

TIDI Choice 917461
White 13" x 18"
2-ply tissue, poly-back,
waffle-embossed bib

U Ultimate HEADREST COVERS

Facial grade tissue provides a pleasing appearance and our softest comfort, along with absorbency and strength. A poly liner also prevents leak-through.



CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
919511	White	10" x 10"	500	FABRICEL® headrest cover, poly-lined
919512	White	13" x 10"	500	FABRICEL headrest cover, poly-lined
919513	White	13" x 13"	500	FABRICEL headrest cover, poly-lined
919514	White	15" x 11-1/2"	500	FABRICEL headrest cover, poly-lined
919515	White	10" x 12", Chayes	500	FABRICEL headrest cover, poly-lined

C Choice HEADREST COVERS

Soft tissue provides comfort for the patient and absorbency for infection control. Poly-lined to prevent leak-through.



CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
919711	White	10" x 10"	500	Tissue/poly headrest cover
919712	Green	10" x 10"	500	Tissue/poly headrest cover
919713	Blue	10" x 10"	500	Tissue/poly headrest cover
919716	Mauve	10" x 10"	500	Tissue/poly headrest cover
919719	Lavender	10" x 10"	500	Tissue/poly headrest cover
919811	White	13" x 10"	500	Tissue/poly headrest cover
919812	Green	13" x 10"	500	Tissue/poly headrest cover
919813	Blue	13" x 10"	500	Tissue/poly headrest cover
919816	Mauve	13" x 10"	500	Tissue/poly headrest cover
919819	Lavender	13" x 10"	500	Tissue/poly headrest cover




E Everyday HEADREST COVERS

A reliable barrier offers peace of mind for patients and protection for equipment.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
919611	White	10" x 10"	500	All-poly headrest cover
919612	Clear	10" x 10"	500	All-poly headrest cover
919621	White	13" x 10"	500	All-poly headrest cover
919622	Clear	13" x 10"	500	All-poly headrest cover
919631	White	13" x 13"	500	All-poly headrest cover

TRAY COVERS

Sized right for your dental tray,
with the right balance of strength
to resist rips and absorbency for
barrier protection.



TIDI Choice
917513 Blue Ritter (B)
tray cover

C Choice

TRAY COVERS

Use our tray covers as a convenient infection prevention barrier. They also minimize noise and improve visibility of tools, burs and files resting on the tray.


CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
917511	White	8-1/2" x 12-1/4"	1000	Ritter (B) tray cover
917512	Green	8-1/2" x 12-1/4"	1000	Ritter (B) tray cover
917513	Blue	8-1/2" x 12-1/4"	1000	Ritter (B) tray cover
917514	Yellow	8-1/2" x 12-1/4"	1000	Ritter (B) tray cover
917516	Mauve	8-1/2" x 12-1/4"	1000	Ritter (B) tray cover
917519	Lavender	8-1/2" x 12-1/4"	1000	Ritter (B) tray cover
917520	TIDI Tooth	8-1/2" x 12-1/4"	1000	Ritter (B) tray cover
917521	White	9-1/2" x 12-3/8"	1000	Chayes (A) tray cover
917523	Blue	9-1/2" x 12-3/8"	1000	Chayes (A) tray cover
917531	White	9-1/2" x 13-1/2"	1000	Midwest (E) tray cover
917533	Blue	9-1/2" x 13-1/2"	1000	Midwest (E) tray cover
917536	Mauve	9-1/2" x 13-1/2"	1000	Midwest (E) tray cover
917541	White	10-1/4" x 15-3/4"	1000	S.S. White (D) tray cover
917543	Blue	10-1/4" x 15-3/4"	1000	S.S. White (D) tray cover
917551	White	11" x 17-1/4"	1000	Weber (C) tray cover
917553	Blue	11" x 17-1/4"	1000	Weber (C) tray cover
917561	White	13-1/2" diameter	1000	Circular (G) tray cover
917571	White	5" x 8"	1000	Mini (F) tray cover
917573	Blue	5" x 8"	1000	Mini (F) tray cover
917581	White	9" x 13-1/2"	1000	Certified tray cover
917591	White	11" x 11"	1000	Cox tray cover

For extra protection and infection prevention, use our tray covers with our poly tray sleeves, found on page 15

Choice COLORS



Choice PATTERNS



TIDI Tooth



SHEATHS + BARRIERS

Aid in the protection from cross contamination for equipment, instruments and other surfaces in the exam environment.



Custom sheath,
Air Techniques
Acclaim camera

CUSTOM INTRA-ORAL CAMERA SHEATHS

Our custom-designed intra-oral camera sheaths fit tightly over the lens for superior image quality and are snug over the body for enhanced productivity. Latex-free and easy to apply and remove.

CODE	QUANTITY	MANUFACTURER	CAMERA MODEL
20917	500 (5 boxes of 100)	Camsight	Optum II
20989	500 (5 boxes of 100)	Camsight	Opticam
20829	500 (5 boxes of 100)	Carestream	STV Pro I & II Kodak 1000
20862	500 (5 boxes of 100)	Cygnus	IDT Oral Vision 5000, 6000
20834	500 (5 boxes of 100)	Dental Vision	Ultracam
20816	500 (5 boxes of 100)	Dentsply	AcuCam 90-degree lens, AcuCam Plus, Polo
20893	500 (5 boxes of 100)	Dentsply	AcuCam Concept III
20856	500 (5 boxes of 100)	Dentsply	Insight
21013	500 (5 boxes of 100)	Digital Doc	Iris
20883	500 (5 boxes of 100)	Digital Doc	NCC Model
20911	500 (5 boxes of 100)	Digital Doc	Aurora
20968	500 (5 boxes of 100)	Digital Doc	Net Design
20990	500 (5 boxes of 100)	Digital Doc	Icon
20810	500 (5 boxes of 100)	DMD	Telicam
20800	500 (5 boxes of 100)	Elca	Elcacam
20818	500 (5 boxes of 100)	Expanded Optics	Explorer
20981	500 (5 boxes of 100)	GC Corp	Crystalcam
20812	500 (5 boxes of 100)	ICS	Vidicam, Intracam, EasyCam, Prospective, Celebrity, Accucure 300
20868	500 (5 boxes of 100)	ICS	Easy Cam, Easy Doc, True Vision
20884	500 (5 boxes of 100)	ICS	Easy Cam II (Tom Cat Supercam and Slimwand, Easy Doc)
20802	500 (5 boxes of 100)	IDT	DT Oral Vision H-3, H-7, H-8, H-9
20987	500 (5 boxes of 100)	Imagin	Hot Rod
21026	500 (5 boxes of 100)	Imagin Systems	Image Master
20991	500 (5 boxes of 100)	Miras	Mirrocam
20984	500 (5 boxes of 100)	Planmeca	AF Camera Sheath
20831	500 (5 boxes of 100)	Pro-Den Systems / Dent-X	Pro-Scope 1000 and Oral Scan, Easy Doc
20806	500 (5 boxes of 100)	Prodentec	Prism
20916	500 (5 boxes of 100)	RF Systems	DP 7S or DP 5/6
20964	500 (5 boxes of 100)	RF Systems	SS21
20867	500 (5 boxes of 100)	Schick	Schick CDR Cam 2000
20855	500 (5 boxes of 100)	Siemens	Ceracam/Minicam Ultra
20970	500 (5 boxes of 100)	Sirona	Sirocam – 3
20899	500 (5 boxes of 100)	Sometech	Dr. Camscope



CUSTOM INTRA-ORAL CAMERA SHEATHS

Continued

CODE	QUANTITY	MANUFACTURER	CAMERA MODEL
20915	500 (5 boxes of 100)	Sometech	Oracam
20928	500 (5 boxes of 100)	Sometech	Oracam
20993	500 (5 boxes of 100)	Sometech	ST-121
20929	500 (5 boxes of 100)	Suni	Suni-Cam II
20969	500 (5 boxes of 100)	Suni	Suni-Cam
20840	500 (5 boxes of 100)	Synca	Logicam / Viper Cam
20909	500 (5 boxes of 100)	Video Dental Concepts	Miniscope & Microscope TD-2000
20963	500 (5 boxes of 100)	Video Dental Concepts	Quickcam Smile
20805	500 (5 boxes of 100)	Vision Tools	Spectra Vu
20852	500 (5 boxes of 100)	Welch Allyn	Reveal
20932	500 (5 boxes of 100)	Win US Tech	WIN 100D
20833	500 (5 boxes of 100)	Yoshida Dental	Pict Cam
20876	500 (5 boxes of 100)	Yoshida Dental	Pict Cam
20905	500 (5 boxes of 100)	Yoshida Dental	Pict Cam

CUSTOM DIGITAL X-RAY SENSOR SHEATHS

Durable, latex-free sheaths have soft edges for patient comfort. A custom fit for your sensor means no excess plastic in the patient's mouth.

CODE	QUANTITY	MANUFACTURER	SIZE
20977	500 (5 boxes of 100)	Carestream/Kodak 6100	Size 0
20978	500 (5 boxes of 100)	Carestream/Kodak 6100	Size 1
20979	500 (5 boxes of 100)	Carestream/Kodak 6100	Size 2
20841	500 (5 boxes of 100)	Cygnus Ray MPS	Size 1
20999	500 (5 boxes of 100)	Dexis	Universal
20957	500 (5 boxes of 100)	Palodex	Size 1
20861	500 (5 boxes of 100)	Regam	Size 2
20998	500 (5 boxes of 100)	Schick	Size 0
20824	500 (5 boxes of 100)	Schick	Size 1
20825	500 (5 boxes of 100)	Schick	Size 2
20819	500 (5 boxes of 100)	Suni	Size 1
20997	Bulk (1 box of 500)	Suni	Extra Large – X-Ray Over Bag
20933	500 (5 boxes of 100)	Progeny	Size 2
20934	500 (5 boxes of 100)	Progeny	Size 1

CUSTOM BARRIERS FROM EQUIPMENT MANUFACTURERS

TIDI custom-fit barriers, designed for a perfect fit, are available direct from the equipment manufacturer. **For barriers designed to fit the following models, contact the manufacturers at the numbers listed below.**

CODE	MANUFACTURER	CAMERA	CONTACT INFORMATION
20927	Air Techniques	Acclaim	Air Techniques Customer Service: 516-433-7676 (Option 1)
20836	Air Techniques	Vistacam	Air Techniques Customer Service: 516-433-7676 (Option 1)
20908	Gendex	AcuCam IV	Imaging Sciences: 800-205-3570
21003	Dexis	DexCam 3	Imaging Sciences: 800-205-3570
21004	Gendex	GXC-300	Imaging Sciences: 800-205-3570
20882	Progeny	OV-1000	Progeny Customer Service: 888-924-3800
20967	Sirona	Sirocam	Sirona USA: 800-659-5977
20892	Sirona	Sirocam 2	Sirona USA: 800-659-5977
20941	Sopro	S717	Sopro Customer Service: 800-289-6367
20972	Sota	Claris	Sota Customer Service: 888-266-7682
20809	Schick Technologies	USB	Schick Customer Service: 1-888-4-SCHICK
21007	Carestream	Kodak 1500	Carestream: 800-944-6365



Soft X-ray sensor material provides comfort for patients.

CUSTOM X-RAY SENSOR SHEATHS

CODE	MODEL	SIZE	CONTACT INFORMATION
20890	Sirona	Size 2	Sirona USA: 800-659-5977
20904	Sirona	Size 1	Sirona USA: 800-659-5977



STANDARD FIT SURFACE BARRIERS

Conveniently protect surfaces against cross contamination.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
915001	Clear	27-1/2" x 24"	225/Box, 4 Boxes/Case	Poly chair cover sleeve, on a roll
915002	Clear	2-1/2" x 10"	500/Box, 24 Boxes/Case	Poly air/water syringe sleeve, closed end
915003	Clear	4" x 6"	500/Box, 24 Boxes/Case	T-style light handle cover
915004	Clear	15" x 26"	250/Box, 4 Boxes/Case	Poly X-ray equipment sleeve
915005	Clear	10-1/2" x 14"	500/Box, 12 Boxes/Case	Poly tray sleeve
915006	Clear	2-1/2" x 10"	500/Box, 24 Boxes/Case	Poly air/water syringe sleeve, open end
915007	Clear	11-5/8" x 14-1/2"	500/Box, 4 Boxes/Case	Poly tray sleeve
915008	Clear	11-1/2" x 16"	500/Box, 4 Boxes/Case	Poly tray sleeve

Simplified surface protection. Additional headrest covers on page 9.



915003 T-style light handle cover



PERSONAL
PROTECTIVE
EYEWEAR



**FLIPEASE®
EYE SHIELD**

- No assembly required
- FLAT for compact storage, FLIP to wear
- Extremely lightweight and comfortable
- Wrap-around design for front and side eye protection
- Optical-grade, anti-static lens

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
9300FE-50	5 assorted colored frames/clear lens	One size fits most	50 (2 boxes of 25)	One piece disposable eyewear; stores flat and assembles with one simple step



Flat, already assembled, and ready-to-wear!



Just flip to wear!




It's that easy!



RESPOSABLES® EYE SHIELD

- Weighs only 3/10 of an ounce!
- Available in pre-assembled, ready-to-wear option or as two separate components (lens and frame)
- “Wrap-around” design provides both front and side eye protection
- Optical-grade replacement lenses; available in clear, UV-tinted and amber-tinted (blue-light absorbent)
- Separately available replacement frames and lenses for maximized cost-effective replenishment

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
9210-50	Clear	One Size	50 (2 envelopes of 25)	Disposable lens for use with 9211-100
9210-100	Clear	One Size	100 (4 envelopes of 25)	Disposable lens for use with 9211-100
9210-250	Clear	One Size	250 (10 envelopes of 25)	Disposable lens for use with 9211-100
9210LT-50	Bronze	One Size	50 (2 envelopes of 25)	Disposable lens for use with 9211-100
9210LT-250	Bronze	One Size	250 (10 envelopes of 25)	Disposable lens for use with 9211-100
9210CL-50	Amber	One Size	50 (2 envelopes of 25)	Disposable lens for use with 9211-100
9210CL-200	Amber	One Size	200 (20 envelopes of 10)	Disposable lens for use with 9211-100
9211-50	Assorted frames	One Size	50 (10 bags of 5 frames)	Reusable frames for use with 9210's
9211-100	Assorted frames	One Size	100 (10 bags of 10 frames)	Reusable frames for use with 9210's
9452	OFFICE PACK	One Size	10 frames and 25 lenses	Disposable eye shields
9210A-50	Assorted color frames/clear lenses	One Size	50 assembled eye shields	Assembled eye shields

**RESPOSABLES® FACE SHIELD**

- Weighs less than 1.5 ounces!
- Separately available replacement visors and lenses for maximized cost-effective replenishment
- Designed to maximize airflow and prevent fogging
- Provides over-the-top, side and front face protection
- Includes a clear, snap-in nosepiece for support
- Disposable, optical grade, anti-static replacement shields

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
2210-50	Clear	13" x 7"	50	Disposable optical grade shield; for use with 2211-100
2210-100	Clear	13" x 7"	100	Disposable optical grade shield; for use with 2211-100
2211-50	Assorted	One size fits most	50	Reusable visor; for use with 2210-100
2211-100	Assorted	One size fits most	100 (10 bags of 10 visors)	Reusable visor; for use with 2210-100
2452	Assorted Visors/ Clear Lenses	One size fits most	3 visors; 15 lenses	Reusable visor and optical grade face shield



Resposables Face Shield-
Patent #6,564,804

DISPOSABLE FACE SHIELD • Preassembled • Adjustable strap for comfortable fit • Latex free

CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
FS850-30	Clear	Standard 11.5" wide x 9" long	30 (1 box of 30)	One piece faceshield with foam headband, elastic strap and optical-grade lens
FS850-120	Clear	Standard 11.5" wide x 9" long	120 (4 boxes of 30)	One piece faceshield with foam headband, elastic strap and optical-grade lens
FS9605-100	Clear	Full Size 14" wide x 11.5" long	100	One piece faceshield with foam headband, elastic strap and optical-grade lens

ACCESSORIES

CODE	COLOR	SIZE	UNIT QUANTITY
9505-50	Neck Strap 5 Assorted Colors	One size	50 Neck Straps
9505-100	Neck Strap Assorted Colors	One size	100 Neck Straps
9810-50	Side Shield Clear	One size	50 (25 pairs)
9902-100	Nose Clips Assorted Colors	One size	100 (50 per bag/2 bags per box) Disposable Nose Clip
9902-1000	Nose Clips Assorted Colors	One size	1000 (50 per bag/20 bags per box) Disposable Nose Clip



FACEMASKS



FACEMASKS

TIDIShield™ fluid resistant facemasks were designed from end user feedback to offer best-available protection and comfort. Manufactured in the USA using soft, breathable materials, our facemasks feature longer, more comfortable earloops and surpass protection standards.

Facemask features:

- Made in USA
- Manufactured of SMS material for greater breathability
- New bright, vibrant color choices
- Longer earloop design for greater comfort and quick/easy donning
- Odorless
- Latex free
- Free of dyes
- Nosepiece is manufactured into facemask folds creating a lighter, softer, more comfortable facemask



CODE	COLOR	SIZE	UNIT QUANTITY	DESCRIPTION
9040DB	Ocean Blue	Universal	10 boxes of 50/ 500 facemasks/case	TIDIShield™ Procedure Facemask
9040DSA	Pink Punch	Universal	10 boxes of 50/ 500 facemasks/case	TIDIShield™ Procedure Facemask
9040DY	Lemon Drop	Universal	10 boxes of 50/ 500 facemasks/case	TIDIShield™ Procedure Facemask

ASTM F2100-07 PERFORMANCE CLASSES FOR FACEMASKS

TIDI facemasks are tested with the ASTM F2100-07 “Standard for Performance of Materials Used in Medical Facemasks.” This is the only specific standard available that provides performance classes for facemasks to fall into. The performance classes range from low, moderate, and high barriers.

REQUIREMENTS	TIDI'S 9040 FACEMASKS
ASTM F1862 Resistance to penetration by synthetic blood (mm HG)	120
ASTM F2101 Bacterial Filtration Efficiency (BFE)	≥ 99
ASTM F2299 Particulate Filtration Efficiency (PFE)	≥ 99
Delta P Breathability	< 2.0
Flame Spread	Class 1

DONNING INSTRUCTIONS



1. Position facemask against face with colored side out.



2. Secure nosepiece comfortably over bridge of nose and pull mask under chin.




3. Pull earloop tight around ears. Confirm mask is secure around circumference of face for customized safety.



CUPS + SALIVA EJECTORS

Our cups' ribbed design offers a sure grip. Flexible saliva ejectors maintain their shape for effective hygiene.



TIDI Choice 917902
Green 18" x 22" 3-ply tissue,
poly-back, contour neck bib.

9220 5-ounce TIDI Tooth
print waxed paper cup.

PATIENT CUPS

Match colors and prints to our tissue and paper products for a friendly and coordinated exam environment.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
9211	White	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9212	Green	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9213	Blue	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9214	Yellow	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9215	Gray	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9216	Mauve	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9217	Beige	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9218	Peach	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9210	Lavender	5 oz.	1000, 100/bag	Plastic cup, ribbed design
9241	White	3.5 oz.	1000, 100/bag	Plastic cup, ribbed design
9243	Blue	3.5 oz.	1000, 100/bag	Plastic cup, ribbed design
9220	TIDI Tooth	5 oz.	1000, 100/bag	Waxed paper cup
9221	TIDI Tooth	5 oz.	500, 100/bag	Waxed paper cup
9226	Jazz	5 oz.	1000, 100/bag	Waxed infused paper cup
9240	Jazz	3 oz.	5000, 100/bag	Waxed infused paper cup
9231	White	6 oz.	2500, 250/bag	Evacuation cup liner, ADEC


COLORS



PATTERNS



Jazz



TIDI Tooth

CUP DISPENSER

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
9246	Chrome	3 or 5 oz.	1	Adjustable cup dispenser

SALIVA EJECTORS

Saliva ejectors contain copper-plated wire, molded in, allowing them to bend easily and then hold their shape. Smooth tips are securely fastened.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
9306	White	6"	1000, 100/bag	White ejector with white tip
9307	Clear	6"	1000, 100/bag	Clear ejector with blue tip
9304	White	Standard	1000, 100/bag	High-volume evacuator tip
9305	White	Standard	2000, 100/bag*	Hygroformic ejector
9333	White	Standard	1000, 100/bag	High-volume evacuator tip, vented

*5 adaptors/100 count bag



ACCESSORIES

COTTON, GAUZE + NONWOVENS

Choose the products that meet your needs, whether pure, natural cotton or nonwovens for maximum absorbency.

COTTON-FILLED GAUZE SPONGES

Cotton-filled sponges feature a woven gauze facing with an absorbent cotton fill, in both sterile and non-sterile packaging.

BRAND	CODE	SIZE	UNIT QUANTITY	DESCRIPTION
VENTURE™	919000	2" x 2"	5000 (25 sleeves of 200)	Gauze sponges with cotton filler, non-sterile
VENTURE	908223	3" x 3"	4000 (40 sleeves of 100)	Gauze sponges with cotton filler, non-sterile
VENTURE	908224	4" x 4"	2000 (20 sleeves of 100)	Gauze sponges with cotton filler, non-sterile
TIDI	919010	2" x 2"	5000 (25 sleeves of 200)	Gauze sponges with cotton filler, sterile



ALL-GAUZE SPONGES

Gauze sponges are ideal for general cleaning and packing.

VENTURE	908282	2" x 2"	5000 (25 sleeves of 200)	8-ply gauze sponges, non-sterile
VENTURE	908294	2" x 2"	8000 (40 sleeves of 200)	12-ply gauze sponges, non-sterile
VENTURE	908206	3" x 3"	4000 (20 sleeves of 200)	8-ply gauze sponges, non-sterile
VENTURE	908293	3" x 3"	4000 (20 sleeves of 200)	12-ply gauze sponges, non-sterile
VENTURE	908292	4" x 4"	4000 (20 sleeves of 200)	8-ply gauze sponges, non-sterile
VENTURE	908295	4" x 4"	2000 (10 sleeves of 200)	12-ply gauze sponges, non-sterile
VENTURE	908291	4" x 4"	2000 (10 sleeves of 200)	16-ply gauze sponges, non-sterile
VENTURE	908200	4" x 8"	2000 (10 sleeves of 200)	12-ply gauze sponges, non-sterile

COTTON ROLLS

TIDI non-linting non x-ray detectable cotton rolls are machine-rolled into soft, pliable cylinders. Venture rolls offer quality at an economical price.

TIDI	969120	1-1/2" x 5/16"	2000 (40 bundles of 50)	Non-sterile, #1 small diameter
TIDI	919121	1-1/2" x 3/8"	2000 (40 bundles of 50)	Non-sterile, #2 medium diameter
TIDI	969122	3/8" x 6"	500 (10 bundles of 50)	Non-sterile, #2 medium diameter
TIDI	969123	1-1/2" x 1/2"	2000 (40 bundles of 50)	Non-sterile, #3 large diameter
VENTURE	919131	1-1/2" x 3/8"	2000 (40 bundles of 50)	Non-sterile, #2 medium diameter
TIDI	969126	1-1/2" x 3/8"	2000 (40 bundles of 50)	Sterile , #2 medium diameter

NONWOVEN SPONGES

Nonwoven sponges are a highly absorbent, low-linting option for prepping and packing. Also convenient for wiping down instruments and workspace.

OMNISORB®	942022	2" x 2"	5000 (25 sleeves of 200)	4-ply (closed weave), non-sterile
OMNISORB	942033	3" x 3"	4000 (20 sleeves of 200)	4-ply (closed weave), non-sterile
OMNISORB	942044	4" x 4"	2000 (10 sleeves of 200)	4-ply (closed weave), non-sterile
VENTURE	912022	2" x 2"	5000 (25 sleeves of 200)	4-ply (closed weave), non-sterile
VENTURE	912033	3" x 3"	4000 (20 sleeves of 200)	4-ply (closed weave), non-sterile
VENTURE	912044	4" x 4"	2000 (10 sleeves of 200)	4-ply (closed weave), non-sterile
VENTURE	260002	2" x 2"	8000 (40 sleeves of 200)	4-ply (open weave), non-sterile
VENTURE	260003	3" x 3"	4000 (20 sleeves of 200)	4-ply (open weave), non-sterile
VENTURE	260004	4" x 4"	2000 (10 sleeves of 200)	4-ply (open weave), non-sterile



DURAWICK™ COUNTER TOWEL

DuraWick™ counter towels were specially created to be used to last through several rounds of tools being placed on them from your ultrasonic cleaning. The top layer of material wicks away moisture, which is then absorbed into the middle layer of the product, keeping your sterilization area clean and sanitary.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
983914	Blue	13" x 18"	100	Wicking towel for ultrasonic tool cleaning, 3-ply Rayon/Tissue/Poly

CSR WRAP

Sterilization wrap conforms easily around various shapes and maintains its wrapped integrity to provide barrier protection against saline, blood and water while allowing effective penetration for a variety of sterilization methods.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
960172	Blue	12" x 12"	1000	Sterilization wrap
960173	Blue	15" x 15"	1000	Sterilization wrap
960174	Blue	18" x 18"	500	Sterilization wrap
960175	Blue	20" x 20"	500	Sterilization wrap
960176	Blue	24" x 24"	500	Sterilization wrap
960177	Blue	30" x 30"	400	Sterilization wrap
960178	Blue	36" x 36"	200	Sterilization wrap

STERILIZATION POUCHES

Self-sealing and have external color indicators for peace of mind.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
5549	Clear one side	2-3/4" x 10"	300	Sterilization pouches, dual indicators, latex free
5550	Clear one side	3-1/4" x 6-1/2"	200	Sterilization pouches, dual indicators, latex free
5551	Clear one side	3-1/2" x 10"	200	Sterilization pouches, dual indicators, latex free
5552	Clear one side	5-1/4" x 11"	200	Sterilization pouches, dual indicators, latex free
5553	Clear one side	7-1/2" x 12"	200	Sterilization pouches, dual indicators, latex free
5543	Blue transparent	2-1/2" x 8-1/4"	1000	Autoclave envelopes, contain latex in adhesive strip
5544	Blue transparent	3-1/2" x 9-3/4"	1000	Autoclave envelopes, contain latex in adhesive strip
5545	Blue transparent	6" x 9-3/4"	1000	Autoclave envelopes, contain latex in adhesive strip
30203	White	2-1/2" x 1-1/2" x 10-1/2"	1000	Autoclave envelopes, contain latex in adhesive strip

AUTOCLAVABLE TOWELS

A poly layer, sandwiched between layers of tissue, allows for use in kits and in autoclaves, while maintaining the integrity of the towels.

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
8253	White	13" x 18"	500	Tissue/poly/tissue autoclavable towel
8252	White	19" x 22"	300	Tissue/poly/tissue autoclavable towel
8251	White	19" x 30"	300	Tissue/poly/tissue autoclavable towel



INFECTIOUS WASTE DISPOSAL

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
8653	Red	24" x 24" (10 gal)	250	High-density (12 microns) infectious waste (Biohazard) bags
8654	Red	24" x 33" (16 gal)	250	High-density (12 microns) infectious waste (Biohazard) bags
8656	Red	33" x 40" (33 gal)	250	High-density (17 microns) infectious waste (Biohazard) bags
8657	Red	40" x 48" (40 gal)	250	High-density (17 microns) infectious waste (Biohazard) bags

CHAIRSIDE DISPOSAL BAGS

CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
950232	White Plain	6-1/2" x 3-1/8" x 11-3/8"	2000	Flame-retardant paper bag, chairside
950231	Blue Floral	6-1/2" x 3-1/8" x 11-3/8"	2000	Flame-retardant paper bag, chairside
950251	Blue Soft Shell	6-1/2" x 3-1/8" x 11-3/8"	2000	Non-flame-retardant paper bag, chairside
950252	White Plain	6-1/2" x 3-1/8" x 11-3/8"	2000	Non-flame-retardant paper bag, chairside
950234	Blue Floral	6-1/2" x 3-1/8" x 11-3/8"	2000 (20 bags of 100)	Non-flame-retardant poly bag, chairside

STORAGE ENVELOPES

CODE	BRAND	SIZE	UNIT QUANTITY	DESCRIPTION
950222	Kraft	8-1/2" x 10-1/2"	1000/case, 20 bundles/case, rotated	X-ray storage envelopes
950221	Kraft	10-1/2" x 12-1/2"	1000/case, 20 bundles/case, rotated	X-ray storage envelopes
950218	Kraft	14-1/2" x 17-1/2"	500/case, 10 bundles/case, rotated	X-ray storage envelopes, heavyweight
950220	Kraft	14-1/2" x 17-1/2"	500/case, 10 bundles/case, rotated	X-ray storage envelopes, lightweight
950219	Kraft	15" x 18"	50/case, bulk packed	X-ray mailer envelopes with insert

KAY-PEES® PROFESSIONAL TOWELS

For sensitive skin, lint-free, low friction.



CODE	COLOR/PRINT	SIZE	UNIT QUANTITY	DESCRIPTION
911898	White	13-1/2" x 17-1/2"	500 (10 boxes of 50)	4-ply tissue, super soft, highly absorbent, edge-embossed
911900	Blue	13-1/2" x 17-1/2"	500 (10 boxes of 50)	4-ply tissue, super soft, highly absorbent, edge-embossed

GLOSSARY

T/P/T: tissue/poly/tissue

A layer of polyethylene film, sandwiched between two sheets of tissue, adds strength and resistance against fluid strike-through to the soft, absorbent tissue.

Tissue (2-ply, 3-ply and 4-ply)

TIDI offers tissue products in configurations to suit a range of needs. Multiple tissue plies add absorbency and strength, while two-ply tissue provides economical softness and absorbency that is well suited for shorter procedures.

Poly

A layer of premium linear polyethylene provides high tensile strength and reliable barrier properties.

FABRICEL®

A proprietary material that fuses a sheet of soft, supple facial-grade tissue and a layer of polyethylene film. This heat lamination process creates a very soft yet super-strong, tear-resistant material that prevents soak-through.

TIDI Custom Barriers

TIDI Products is proud to partner with manufacturers of dental and medical devices that require a custom-fit infection control barrier. Our design team works with manufacturers to create custom sheaths that aid in the protection from cross contamination while allowing healthcare professionals to use the instrument without the barrier interfering with the procedure. If you are a manufacturer and would like more information on how we can develop a custom barrier for your device, please **contact TIDI Customer Service at 800-521-1314.**



DENTAL TRADE ALLIANCE



TIDI supports the Forest Stewardship Council and responsible forestry.

Certified in ISO 9001 and EN/ISO 13485



TIDI Products, LLC

570 Enterprise Drive
Neenah, WI 54956 USA

Customer Service

800.521.1314

+1 920.751.4138

800.837.7770 FAX

www.tidiproductions.com

